2022-10
802020
7.0
M

 TRACOE medical GmbH
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200
info@tracoe.com, www.tracoe.com

TRACOE®
kids

CE
0 1 2 3

**TRACOE® silcosoft® Tracheostomy Tubes for
Neonates / Infants and Children**

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Image 4: Plastic obturator


Image 5: Metal obturator


Image 1: Uncuffed model with plastic obturator



Image 2: Cuffed Model with proximate extension and metal obturator

REF 360, 362

Size	ID mm	OD mm	C mm	TL mm	OL mm	OB mm	BA
2.5	2.5	4.4	26.8	30	55.5	32	120°
3.0	3.0	4.9	28.0	32	57.5	34	120°
3.5	3.5	5.4	29.0	34	59.5	37	120°
4.0	4.0	6.0	29.0	36	61.5	39	120°

REF 361, 363

Size	ID mm	OD mm	C mm	TL mm	OL mm	OB mm	BA	PL
2.5	2.5	4.4	26.8	30	80	32	120°	15
3.0	3.0	4.9	28.0	32	84	34	120°	17
3.5	3.5	5.4	29.0	34	89	37	120°	20
4.0	4.0	6.0	29.0	36	93	39	120°	22


REF 370, 372

Size	ID mm	OD mm	C mm	TL mm	OL mm	OB mm	BA
2.5	2.5	4.4	29.0	38	63.5	40	120°
3.0	3.0	4.9	29.0	39	64.5	41	120°
3.5	3.5	5.4	29.0	40	65.5	43	120°
4.0	4.0	6.0	29.0	41	66.5	44	120°
4.5	4.5	6.6	29.0	42	67.5	45	120°
5.0	5.0	7.3	29.0	44	69.5	48	120°
5.5	5.5	8.4	29.0	46	71.5	50	120°


REF 371, 373

Size	ID mm	OD mm	C mm	TL mm	OL mm	OB mm	BA	PL
2.5	2.5	4.4	29.0	38	93	40	120°	20
3.0	3.0	4.9	29.0	39	94	41	120°	20
3.5	3.5	5.4	29.0	40	100	43	120°	25
4.0	4.0	6.0	29.0	41	101	44	120°	25
4.5	4.5	6.6	29.0	42	107	45	120°	30
5.0	5.0	7.3	29.0	44	109	48	120°	30
5.5	5.5	8.4	29.0	46	111	50	120°	30

REF 360, 362, 370, 372

REF 361, 363, 371, 373

ID: inside diameter (clear width) of tube; **OD:** outside diameter of tube;
C: length of bend; **TL:** length along center line from start of neck flange to bottom of tube; **OB:** length along outer bend from start of neck flange to bottom of tube; **BA:** bending angle; **PL:** proximal longer; **OL:** overall length along center line from bottom of tube to connector end.

DE / Gebrauchsanweisung beachten • **EN** / Please read Instructions for Use • **FR** / Respecter le mode d'emploi • **IT** / Consultare le istruzioni per l'uso • **ES** / Síga las iZstrucciones de uso • **PT** / Ter atenção às instruções de utilização • **DA** / Følg brugsanvisningen • **FI** / Noudata käyttöohjetta • **NO** / Følg bruksanvisningen • **HU** / Olvassa el a használati útmutatót • **SV** / Se bruksanvisningen • **NL** / Rekening houden met de gebruiksaanwijzing • **HR** / Pridržavajte se uputa za uporabu • **SL** / Upoštevajte navodila za uporabo • **CS** / Respektujte návod k použití • **RO** / Respectați instrucțiunile de utilizare • **TR** / Kullanma talimatına dikkat ediniz • **RU** / Соблюдать инструкцию по применению • **PL** / Zastrzyj do instrukcji używania • **EL** / Ακολουθήστε τις οδηγίες χρήσης • **KO** / 사용 설명서를 준수하십시오 • **ZH** / 注意使用说明书 • **AR** / لامل عتسالا ليلد ةاعارم بجي

DE / Achtung. Gebrauchsanweisung beachten • **EN** / Warning. Please read Instructions for Use • **FR** / Attention. Respecter le mode d'emploi • **IT** / Attenzione. Consultare le istruzioni per l'uso • **ES** / Atención. Síga las instrucciones de uso • **PT** / Atenção. Ter atenção às instruções de utilização • **DA** / Vigtigt: Følg brugsanvisningen • **FI** / Huomio. Noudata käyttöohjetta • **NO** / OBS. Følg bruksanvisningen • **HU** / Figyelem. Olvassa el a használati útmutatót • **SV** / OBS! Se bruksanvisningen • **NL** / Let op. Rekening houden met de gebruiksaanwijzing • **HR** / Upozorenje. Pridržavajte se uputa za uporabu • **SL** / Pozor. Upoštevajte navodila za uporabo • **CS** / Pozor. Respektujte návod k použití • **RO** / Atenție. Respectați instrucțiunile de utilizare • **TR** / Dikkat. Kullanma talimatına dikkat ediniz • **RU** / Внимание! Соблюдать инструкцию по применению • **PL** / Ostrzeżenie. Zastrzyj do instrukcji używania • **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης • **KO** / 주의. 사용 설명서를 준수하십시오 • **ZH** / 注意。注意使用说明书 • **AR** / ةاعارم بجي .هيبنت لامل عتسالا ليلد

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Instructions for Use for

TRACOE® silcosoft® Tracheostomy Tubes for Neonates / Infants and Children

Note:

- Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own and your patients' safety, please observe the following safety information.

- The illustrations to which the text refers can be found on the (foldout) illustrated pages at the beginning of these instructions. The numbers and letters in brackets refer to the respective illustrations and product components of the tracheostomy tube. The symbols and icons used are explained on pages 1 – 5.

1. Indications for Use / Intended Use

The TRACOE silcosoft tracheostomy tube is intended to provide direct airway access for a tracheotomized patient up to 29 days. It may be reprocessed for single patient use up to 7 times.

Patient population: Neonates, infants, children, and adolescents.

2. General Description

TRACOE silcosoft is a silicone tracheostomy tube which provides an artificial airway to the lower respiratory tract for neonates, infants, children, and adolescents.

The TRACOE silcosoft tracheostomy tube is provided in a sterile blister pack together with an obturator, a fabric neck strap and a disconnect wedge. The silicone tracheostomy tube is radiopaque, and wire reinforced. It is available in different diameters and lengths, and in cuffed or uncuffed models. The cuffed model is provided with the cuff deflated and flush with the tube. The obturator is made of plastic or metal depending on the type of tube (please see Table 1 on the fold-out pages). The appropriate diameter and length of the tube is to be determined by the physician.

For insertion, the neck is extended with the tracheostoma open and free of obstruction. The tracheostomy tube (with the obturator inside), is inserted directly into the tracheostoma until the neck plate is in contact with the skin surface. Once in place, the obturator is immediately removed and when applicable, the cuff is filled with sterile water. With the tracheostomy tube in position it is secured in place with a neck strap.

The TRACOE silcosoft is intended for single patient use up to 29 days and may be reprocessed 7 times within 29 days. The frequency of reprocessing is based on the physician's care instructions. The device is applicable for mechanically ventilated or spontaneously breathing

patients in hospitals, pre-hospital (EMS), extended care facilities, outpatient clinics, or home care and can be used by individuals trained in tracheostomy care.

This product is supplied with an information card, which includes two detachable labels with product specific details and MR imaging information. The labels facilitate reordering and can be attached to the patient record.

The Images 1 and 2 on the foldout pages represent TRACOE silcosoft models (description see Table 2).

1a	Plastic obturator	5b	Check valve
1b	Metal obturator	5c	Female Luer connector
2	Tube (distal end)	6	15 mm connector
3	Neck flange (plate)	7	Tube (proximal extension)
4	H ₂ O Cuff	7a	Connector hub (silicone)
5	Inflation (filling) line	8	Disconnect wedge
5a	Pilot balloon	9	Neck strap

An overview of the different types and their components is provided in Table 1 also on the foldout pages.

Obturator:

The obturator (**1**) has a handle on the proximal end and a rounded, conically shaped tip on its distal end. The obturator shape is adapted to the tube's inner diameter, length, and bending angle.

- The tracheostomy tube without proximal extension has a plastic obturator (**1a**) with a ring handle, to assist with insertion. For correct alignment of the obturator the display for REF-No on the ring handle should run horizontal, see Image 4.
- The tracheostomy tube with a proximal extension has a metal obturator with a flat handle, to assist with insertion. For correct alignment of the obturator the arrow on the handle shows upwards and the writing runs horizontal, see Image 5.

Tracheostomy Tube:

- The distal end (inside the patient) of the tracheostomy tube (**2**) is wire reinforced, curved and features a smooth round tip.
- The tracheostomy tube includes a fixed neck flange (**3**) which limits the insertion depth and includes eyelet holes for attachment of the neck strap (**9**).
- Tracheostomy tubes without proximal extension have a silicone hub (**7a**) in combination with a neck flange (**3**). The 15 mm connector (**6**) is directly attached to the silicone hub (**7a**).
- Tracheostomy tubes with proximal extension (**7**) have a wire-reinforced tube between neck flange (**3**) and silicone hub (**7a**). The silicone hub provides a smooth, stable transition between the

tube and the 15 mm connector.

- The 15 mm connector (6) provides a smooth surface for connecting the tracheostomy tube to external devices e.g. for mechanical ventilation.

H₂O Cuff Models:

- The silicone cuff (4) is inflated with sterile water, and located on the distal end of the tracheostomy tube and. It is directly connected to the inflation line (5), which is incorporated into the neck flange (3).
- The proximal end of the inflation line includes a pilot balloon (5a), an incorporated self-sealing check valve (5b), with a female Luer connector (5c).

Disconnect Wedge:

- The plastic disconnect wedge (8) is a tool to disconnect facilitate disconnection of external devices from the 15 mm connector of the tracheostomy tube.
- The disconnect wedge slides between the female 15 mm connector of the external device and the male 15 mm connector of the tracheostomy tube to easily separate the devices, see Image 3.

Neck Strap:

- The neck strap (9) is a soft strip of padded fabric that wraps around the patient's neck.
- The ends of the strap include hook and loop fasteners that are inserted through the eyelet holes of the neck flange (3) to secure the tracheostomy tube position.

3. MRI Safety Information

Non-clinical testing has demonstrated the TRACOE silcosoft tracheostomy tube is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla (T) or 3.0 Tesla (T).
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg and a maximum whole head specific absorption rate (SAR) of 3.2 W/kg
- Quadrature driven transmit body coil only.
- The neck flange (3) must be secured in place with the neck strap (9).
- For tubes with H₂O Cuff only: The check valve (5b) must be secured to the skin with medical tape, away from the area of MR diagnostic interest.

In non-clinical testing, the image artifact caused by the tube extends (radially) up to 19 mm from the TRACOE silcosoft tracheostomy tube (non-cuffed) when imaged with a gradient echo pulse sequence and a 1.5T MR system, and up to 28 mm when imaged with a spin echo pulse sequence in a 3.0T MR system.

Tubes with H₂O Cuff: In non-clinical testing, the image artifact, caused by the check valve, extends (radially) up to 107 mm from the check valve when imaged with a gradient echo pulse sequence and a 1.5T MR system, and up to 113 mm when imaged with a spin echo pulse sequence in a 3.0T MR system. Therefore, when MR imaging a tube with the H₂O Cuff, it is recommended to tape the check valve to the patient's skin away from the area of interest.

Warning:

When used in MR Imaging:

- Securely fasten the tube, with the neck strap, to prevent possible movement while in the MR environment.
- Securely affix the check valve (**5b**) away from the area of interest with standard medical tape to prevent movement within the MR environment.
- MR image quality may be compromised if the area of interest is close to the position of the tube and/or inflation valve.

4. Contraindications

- The TRACOE silcosoft tracheostomy tube cannot be used in conjunction with heat emitting devices, e.g. lasers to prevent the tube material from being damaged.
- The H₂O Cuff should not be inflated, when a speaking valve is used and vice versa.

5. General Precautions

- A spare tracheostomy tube should remain at the patient's bedside and available for replacement, as required. This device should be clean and dry.
- Ensure that an alternative airway is always available.
- The disconnect wedge (supplied with the device) should remain at the patient's bedside and be available to assist in the disconnection of external devices, e.g. mechanical ventilation.
- The sterile packaging should be inspected for damage prior to opening. If the packaging is damaged the device should not be used.
- The tracheostomy tube should be inspected for tube integrity and function prior to use/insertion. Verify that the tube is free of obstruction, H₂O Cuff material is not brittle or torn and can be inflated/deflated, absence of kinks, tears or cuts, stable connection between the tube and the neck flange, etc. If the product is damaged, it should be discarded and replaced with a new product.
- Do not use unnecessary force on the tracheostomy tube when connecting to, or disconnecting from, external devices. This may result in damage to the tracheostomy tube, decannulation and/or increased air resistance. Always hold the tracheostomy tube at the base of the 15 mm connector (**6**) when disconnecting from external devices. If external devices are tightly fixed, use the disconnect wedge (**9**) to separate them from the 15 mm connector tube.

- When the tracheostomy tube is used together with other medical devices, please acknowledge their instructions for use. Contact the manufacturer if there are any questions or assistance required.
- The use of tracheostomy tubes can lead to pressure points and skin irritation (e.g. due to moisture). To avoid this place a dressing underneath the neck flange.
- Ensure that the syringe for filling the H₂O Cuff is clean and free of contaminants. The syringe must be removed immediately after use. Failure to do so may result in a deflation of the H₂O Cuff filling system which will reduce efficient ventilation or protection from aspiration.
- To avoid damage to the H₂O Cuff and improve ease of insertion, always ensure that the H₂O Cuff is completely deflated prior to insertion.
- To avoid damage to the H₂O Cuff material it should not be in contact with lidocaine-containing aerosols or ointments.
- When repositioning the patient, ensure that the patient does not squeeze the pilot balloon (5a). This can increase the cuff pressure and therefore harm the trachea.
- Prolonged storage of the tube after reprocessing can increase the risk of bacterial contamination.
- For appropriate handling of the tracheostomy tube follow the standard guidelines for prevention of infection as stipulated by the Center for Disease Control and Prevention (USA), the Robert Koch Institut (DE), or the local or national authority.

6. Warnings

- Do not use, if the sterile packaging has been compromised/damaged, e.g. open edges, holes in packaging etc.
- Improper storage conditions may result in tube damage or contamination.
- Use only water-soluble lubricating gel.
- To prevent an increase in air resistance, ensure that the tube does not become obstructed when applying lubricating gel to the obturator tip.
- Check the position and function of the tube following insertion. Incorrect placement may result in permanent damage to the tracheal mucosa or increased air resistance.
- Do not move or shift the tube once it is in position, as this may damage the stoma / trachea or lead to insufficient ventilation.
- A subtle increase in air resistance is possible, if the tube is not suctioned regularly.
- Use only suction catheters to clear the secretions from the patient's respiratory tract and the tracheostomy tube. Instruments with a conical tip may wedge in the tube and restrict ventilation.
- If the suction catheter gets stuck in the tracheostomy tube rotate the catheter.
- Regularly check that all connections are secure to prevent an inadvertent disconnection of the tube from external equipment and ensure efficient ventilation.

- Only fill the H₂O Cuff with sterile water. Do not fill the H₂O Cuff with air as this would impair the sealing function of the cuff which may result in aspiration or VAP (ventilator associated pneumonia).
- The H₂O Cuff is an HPLV (high pressure low volume) cuff therefore, the measured cuff pressure is not related to the actual cuff pressure exerted on the tracheal wall.
- Before deflating the H₂O Cuff (removal of sterile water), ensure that the patient's upper respiratory tract is unobstructed. When applicable, clear the upper respiratory tract of any secretions through suction and if possible patient coughing.
- Do not overfill the H₂O Cuff, as this poses a risk of permanent damage to the trachea.
- Insufficient filling of the H₂O Cuff could result in an increased risk of aspiration or VAP (ventilator associated pneumonia).
- If the cuff leaks during use, replace the tube and contact TRACOE Customer Service.
- To prevent damage to the stoma or trachea, ensure that the H₂O Cuff is completely deflated (empty) prior to insertion or removal of the tube (see chapter 8). If it is not possible to deflate the cuff, cut the inflation line with a pair of scissors and drain the liquid. With this event, the product is defective and must be replaced.
- Additional proximal extension pieces, i.e. TRACOE *kids* extension piece REF 356 are not allowed to be used in combination with an proximal extended TRACOE silcosoft tube (REF 361, REF 363, REF 371, REF 373).

7. Adverse Reactions

Typical adverse effects during usage of tracheostomy tubes include pressure points, stenosis, skin irritation, a need to cough, and swallowing difficulties.

Irritation, a need to cough or bleeding may also occur during insertion and removal of the tube. In the event of an adverse effect, please contact a medical professional immediately.

8. Functional Description

8.1 Preparing the Tube

Note: If the tube fails any step of inspection, repeat the procedure with a new TRACOE silcosoft tracheostomy tube. Keep the tube with the batch number when contacting TRACOE Customer Service.

The following functions must be checked immediately prior to use:

- Inspect the sterile packaging to ensure it is secure, undamaged and all components are present.
- Open the package and inspect the device for damage prior to use.
- Ensure the obturator can be easily moved in and out of the tube.
- Ensure the envisaged suction catheter is of appropriate size and can be easily inserted through the tube.
- If obturator or suction catheter can not be moved easily, please

ad some water-soluble lubricant gel to the tip of the obturator or catheter.

If the device fails the initial inspection, discard and repeat the procedure with a new device.

Additional inspection requirements for tubes with H₂O Cuff: Inspect the pilot balloon prior to use and perform a H₂O Cuff-filling test to ensure that the cuff is leak-proof and fills uniformly.

H₂O Cuff-filling test: Use a sterile syringe (0 – 10 ml) with ml scale. Fill the syringe with sterile water and connect it to the female Luer connector (5c) of the filling tube (check valve). Press the plunger of the syringe and fill the H₂O Cuff with 2 ml (maximum) of liquid. Inspect the balloon to ensure it is leak proof and fills uniformly.

When the inspection is complete, remove the sterile water from the cuff, by pulling back the syringe plunger. Once all the water is removed the tube (with deflated cuff) is ready for insertion, see chapter 8.4.

Warning:

- If the H₂O Cuff leaks, do not use the tube.
- If the H₂O Cuff cannot be filled or adheres (“sticks”) to the tracheostomy tube, carefully press the cuff between thumb and forefinger to separate it from the tube, then repeat the filling procedure. If this procedure has no effect, do not use the tube.

8.2 Preparing the patient

Ensure that the patient is optimally pre-oxygenated before insertion or re-insertion. To facilitate insertion, slightly overextend the patient's neck, if possible.

Should complications arise during insertion or re-insertion, appropriate safety measures should be taken (e.g. keep on hand a tracheal dilator or a tube with a smaller diameter).

8.3 Inserting the tube

Prior to insertion, verify that the obturator tip has not been inadvertently pushed inside the distal end of the tube. This can be avoided by simultaneously holding the tube and obturator with one hand.

With the obturator properly inserted into the tube, apply a pea-sized amount of water-soluble lubricating gel on the tip without blocking the tube opening.

Extend the patient's neck, and with the stoma clear, insert the tube and obturator (with lubricated tip), into the patient's stoma until the neck flange is in contact with the skin surface. Immediately pull the obturator out of the tube to open the airway through the tube.

8.4 Following Tube Insertion

Secure the tube with the neck strap by attaching one end of the neck strap (9) through the opening on the neck flange (3) and wrapping the strap around the patient's neck, attaching the opposite end of the neck strap through the opposite neck flange opening. Ensure that the strap is

tight, and the tube is secured to prevent movement during clinical use.

For tracheostomy tubes with H₂O Cuff:

- When the tube is secured, fill the H₂O Cuff with sterile water through the pilot valve, using a syringe.
- Fill the syringe with 5 ml of sterile water.
- Use either the minimal leak technique or minimal occlusive volume technique when filling the cuff, to ensure appropriate tracheal seal.
- Record the individual patient-related volume used to fill the H₂O Cuff.
- Verify the position and function (see chapter 6 Warnings) of the tube. Do not adjust the tube position when the H₂O Cuff is filled, as this may damage the stoma and/or tracheal mucosa.

8.5 Connecting/Disconnecting External Equipment

To connect to external equipment (e.g. mechanical respirator): Firmly hold the base of the 15 mm connector and gently twist the connection end of the external device until it is securely attached to the tracheostomy tube. If in doubt, twist the connection end on and off several times, in order to confirm the amount of force needed to ensure the connection is secure and the external device can be easily disconnected at a later time.

When disconnection of the external device is difficult, use the disconnect wedge. Slide the opening of the disconnect wedge between the 15 mm connector and external device until the two devices separate (see Image 3).

8.6 Deflating the H₂O Cuff

Before deflating the cuff, ensure that secretions cannot enter the lower respiratory tract.

Attach a syringe (with the plunger pushed in) to the female Luer connector (5c). When attached, pull the plunger back until all liquid is removed from the cuff. The H₂O Cuff must be completely deflated (empty) prior to removal.

8.7 Removing the Tube

- Disconnect the neck strap from the neck flange.
- Slightly overextend the patient's neck, if possible.
- For tubes with H₂O Cuff, ensure the cuff is completely empty of liquid.
- Firmly hold the base of the 15 mm connector and gently pull the tracheostomy tube from the stoma.
- Following removal, the tube should be cleaned as soon as possible to prevent encrustation of fluids.
- If the product is damaged, do not reuse the tube and inform TRACOE Customer Service.

9. Cleaning, Reprocessing and Storage

Following the cleaning or reprocessing procedure, a visual inspection prior to re-insertion should be performed, see chapter 8.1.

CAUTION:

- Cleaning and reprocessing the device outside these guidelines, or failure to clean the device properly can result in damage to the tube, possible obstruction or increased air resistance and damage to the tracheal mucosa, infection or irritation/inflammation of the tracheal stoma.
- Reprocessing / sterilization methods other than the ones described in section 9.2 are not permitted, such as deep, pulse or intermittent vacuum sterilization, EO, H₂O₂ and O₃, microwave sterilization, or Dry Heat.
- To avoid recontamination after cleaning and disinfection, only hold the tube at the base of the 15 mm connector and not at the insertion end of the tube

9.1 Cleaning

Cleaning of the device is intended to remove any bodily fluids or encrustation that may inhibit its clinical use.

The frequency of cleaning must be defined by the physician.

The tube (2) and obturator (1) must be cleaned before reprocessing (see section 9.2), and/or reuse.

The following instruction for manual cleaning applies to all TRACOE silcosoft models and sizes:

The tracheostomy tube should be cleaned immediately after removal from the stoma to prevent drying of soil and contaminants.

To clean the tube and obturator, rinse the devices separately under lukewarm (40 °C/104 °F) potable water at least for 3 minutes, until they are visibly clean and free of contaminants/encrustations. Particular attention should be taken to ensure the inside of the tube is thoroughly rinsed. If encrustations are still present after rinsing, soak the device in warm water (up to 65 °C/149 °F) up to 30 minutes and rinse again. Repeat these steps as required until the tube is visually clean. All areas of the tube and obturator should be inspected, in adequate light, to ensure the device is free of contaminants and encrustations. When cleaning, take care not to damage the H₂O Cuff.

Following the cleaning process, place the tube and obturator on a clean dry towel and air-dry in an area free of airborne contaminants. The tube and obturator are considered dry when there is no visual evidence of residual water.

Following the drying process, the tube can be reinserted, or sterilized with steam sterilization, see 9.2.2.

If, after cleaning, the tube and obturator will be disinfected with hot

water disinfection (9.2.1), drying is not required.

CAUTION:

Do not use any cleaning agents, brushes or swabs during cleaning as they can damage the tube.

9.2 Reprocessing

The TRACOE silcosoft is intended for single patient use up to 29 days and may be reprocessed 7 times within 29 days. The frequency of reprocessing is based on the physician's care instructions.

NOTE:

The disconnect wedge and neck strap cannot be reprocessed.

9.2.1 Hot Water Disinfection (typically used in home care)

Fill a pot (with lid) with potable water and place over a heat source until the water is boiling (100 °C/212 °F). Place the clean tube and obturator, separately, into the pot of boiling water. Immediately cover the pot with a lid and allow the water to continue boiling for 3 minutes. Then carefully remove the pot from the heat source. The tube and obturator should remain in the covered pot of water for 30 minutes, at which time the water should be cool. Remove the tube by holding the 15 mm connector end (6), and the obturator by the handle (1).

Place the tube and obturator on a clean dry towel and air-dry in an area free of airborne contaminants. The tube and obturator are considered dry when there is no visual evidence of residual water.

All areas of tube and obturator should be visually inspected, in adequate light, to ensure the device is free of contaminants and there is no visible damage to the devices, e.g. cracks, holes etc. Following inspection, the obturator can be inserted in the tube and the device reinserted into the patient following the instruction defined in section 8.3.

If the device is not immediately reinserted into the patient, place in a clean dry container that can be sealed (e.g. lid) to reduce the risk of contamination.

CAUTION:

- Risk of scalding with hot water.
- Make sure the tube is not hot when re-inserted into the stoma.

9.2.2 Autoclave Sterilization (Moist Heat/Steam)

Place the clean, dry tube and obturator in separate sterilization pouches (use only an approved sterilization pouch).

Place the tube and obturator for 30 minutes at 121°C (250 °F) in a gravity displacement steam sterilizer, with a minimum drying time of 15 minutes. Care must be taken to follow the steam sterilizer's manufacturer instructions for use. Once the device has been sterilized, it can be re-inserted into the patient as defined in section 8.3.

If the sterilized tube is not used immediately, store the sterilized tube and obturator in their sterilization pouch, as defined in 9.3.

CAUTION:

- The device is single patient use; therefore the tube must be returned to the same patient.

9.3 Storage

- a) Store the factory delivered tracheostomy tubes in their original packaging in a dry location, away from sunlight until initial use.
- b) Store reprocessed tracheostomy tubes in a clean covered container or sterilization pouch, in a clean dry location, away from sunlight. Re-insert as soon as possible.
- c) Improper storage conditions may result in tube damage or contamination. Do not store the reprocessed devices for more than 29 days or beyond the time limit defined by the sterilization pouch.

10. Packaging

TRACOE silcosoft is provided sterile (sterilized with ethylene oxide) which allows the product to be inserted under sterile conditions, e.g. initial insertion in surgical suite (operating room). TRACOE silcosoft does not require a sterile environment during normal use or reprocessing. Sterility is only ensured if the sterile device packaging is intact and sealed on all sides prior to initial use.

11. Disposal

Used products are to be disposed of in accordance with the valid national regulations, or clinical procedures governing biohazardous waste materials. Ensure the device has been rinsed and free of contaminant before disposal, see section 9.1.

12. Returns

Returned products, that have been used, will only be accepted subject to prior consultation and if a completed decontamination certificate and complaint report is enclosed. These forms are available either directly from TRACOE medical or via the website www.tracoe.com.

13. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be affected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on our website at www.tracoe.com.

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